

Human CD40 ligand recombinant protein

PX-P4016-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	17.95kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	0,8499999999999999
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	50 mM Tris-HCl, pH 8.0, 150 mM NaCl 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MENSFEMQKGDQNPQIAAHVISEASSKTTSVLQWAEKGYTMSNNLVTLENGKQLTVKRQG
LYYIYAQVTFCNSNREASSQ
APFIASLCLKSPGRFERILLRAANTHSSAKPCGQQSIHLGGVFELQPGASVFNVTDPQSQVSHGT
GFTSFGLLKLLEHHH
HHH

More info

UNIPROT ID

P63304

UNIPROT LINK

<http://www.uniprot.org/uni-prot/P63304>

SYNONYMS

CD40L, CD154, TRAP,
HIGM1, IGM, IMD3,
TBAM, T-BAM, TNFSF5,
Gp39, TNF Superfamily
Member 5, Hyper-IgM
Syndrome, TNF-Related
Activation Protein, T-Cell
B-Cell Activating Molecule

For research use only.