

Human CD152/Cytotoxic T-lymphocyte protein 4 recombinant protein

PX-P4018-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	18.51kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,94999999999999996
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	50 mM Tris-HCl, pH 8.0, 150 mM NaCl 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MACLGFQRHKAQLNLA TRTPCTLLFFLLFIPVFKAMHVAQPAVV LASSRGIASFVCEYASP
GKATEVRVTVLRQADSQVTEVCAATYMMGNELTFLDDSICTGTSSGNQVNLTIQGLRAMDTG
LYICKVELMYPPPYLIGINGTQIYVIDPEPCPDSGDGSHHHHHH

More info

UNIPROT ID

P16410

UNIPROT LINK

<http://www.uniprot.org/uni-prot/P16410>

SYNONYMS

CTLA4, Cynw2, CD152, CELIAC3, CTLA4, GSE, IDDM12, Celiac Disease 3, Cytotoxic T-lymphocyte protein 4