

## Human IL33 isoform 1 partial (112-270) recombinant protein

PX-P4020-10

### OVERVIEW

---

<b>SIZE</b>	10 ug
<b>ORIGIN SPECIES</b>	Human
<b>FRAGMENT</b>	
<b>PROTEIN DELIVERED WITH TAG</b>	Yes
<b>MOLECULAR WEIGHT WITH TAG IF ANY</b>	19.92kDa
<b>DELIVERY CONDITION</b>	any condition

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	E.coli
<b>PURITY</b>	0,8499999999999998
<b>PROTEIN ACCESSION</b>	
<b>FORM</b>	Lyophilized
<b>BUFFER</b>	PBS, pH 7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

### PROTEIN SEQUENCE

---

MGSHHHHHHSGDDDDKMSITGISPITEYLASLSTYNDQSITFALEDESYEIYVEDLKKDEKKDK  
VLLSYYESQHPSNESGDGVDGKMLMVTLSPTKDFWLHANNKEHSVELHKCEKPLPDQAFFVL  
HNMHSNCVSFECKTDPGVFIGVKDNHLALIKVDSSSENLCCTENILFKLSET

### More info

#### UNIPROT ID

O95760

#### UNIPROT LINK

<http://www.uniprot.org/uni-prot/O95760>

#### SYNONYMS

DV27, C9ORF26, IL1F11, NFHEV, Interleukin-1 Family, Member 11, Nuclear factor from high endothelial venules, Interleukin 33, IL33,IL33 isoform 1 partial (112-270)