

## Recombinant Human C reactive protein (CRP)

PX-P4056-10

### OVERVIEW

---

<b>SIZE</b>	10 ug
<b>ORIGIN SPECIES</b>	Human
<b>FRAGMENT</b>	
<b>PROTEIN DELIVERED WITH TAG</b>	Yes
<b>MOLECULAR WEIGHT WITH TAG IF ANY</b>	unknown
<b>DELIVERY CONDITION</b>	any condition

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	mammalian
<b>PURITY</b>	0,90000000000000002
<b>PROTEIN ACCESSION</b>	
<b>FORM</b>	Lyophilized
<b>BUFFER</b>	50 mM Tris-HCl pH 7.5, 150 mM NaCl 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>STABILITY &amp; STORAGE</b>	

### PROTEIN SEQUENCE

---

/

### More info

#### UNIPROT ID

P02741

#### UNIPROT LINK

<http://www.uniprot.org/uni-prot/P02741>

#### SYNONYMS

CRP, C-RP, PTX1,  
Pentraxin-Related,  
Pentraxin 1