

Human TNFSF9/4-1BBL recombinant protein

PX-P4060-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	22.64kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,5
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

MGSHHHHHHSGACPWAVSGARASPGSAASPRLLREGPELSPDDPAGLLDLRQGMFAQLVAQN
VLLIDGPLSWYSDPGLAGVSLTGGLSYKEDTKELVVAKAGVYYVFFQLELRRV VAGEGSGSV
SLALHLQPLRSAAGAAALALTVDLPPASSEARNSAFGFQGRLLHLSAGQRLGVHLHTEARARH
AWQLTQGATVLGLFRVTPEIPAGLPSRSE

More info

UNIPROT ID

P41273

UNIPROT LINK

<http://www.uniprot.org/uni-prot/P41273>

SYNONYMS

4 1BB L, 4 1BB ligand, 4 1BBL, 4-1BB ligand, 4-1BBL, Cd137L, Cd157L, Homolog of mouse 4 1BB L, Homolog of mouse 4 1BBL, ILA ligand (TNF related), Ly63L, Receptor 4 1BB ligand, TNF superfamily member 9, TNFL9_HUMAN, Tnfsf9, TNLG5A, Tumor necrosis factor (ligand) superfamily member 9, Tumor necrosis factor ligand 5A, Tumor necrosis factor ligand superfamily member 9, Tumor necrosis factor superfamily member 9