

Bifidobacterium longum Amidase Recombinant Protein

PX-P4064-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Bifidobacterium longum
FRAGMENT	Partial
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	46.76kDa
DELIVERY CONDITION	

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	0,6999999999999996
PROTEIN ACCESSION	
FORM	
BUFFER	
STABILITY & STORAGE	

PROTEIN SEQUENCE

MGCVGGGLAFSSAQAQADTYSDLINAQNQHAASVQREAEKQQLAGASQDLANKVLELDDL
TNNKIVAAQAKVTQANEDAATAQDEADAASGRLSAAQKDKETLEEQIKQTGKDYDDAHAAV
AQLARDEMHGNSASDVMSVVTGATSTQDFVNSMQSRDALSRNEANAASSAATSLSTSKNRG
ERLAAIEKQIAVLKTADEKAASAQTAAETAQSERDALDKLRQEGEARRELSMIDSLDSQS
AKQAAQTVLIASQVDSYNRQFQKEQQDAANRVDTGNQGGTPSTPVPAPAPAPAPAPAPAPA
PAPSVGGQGTNGDYGWAYATGQCTYWAYERRRQMIGITPSYLGNGGDWWRNAPSGLRV
DHNQVGAALSFLPGQDGADGTWGHVAVVEAVYGDGTFQISEMNVGGLWMMNYRTLNL
GQYWFVHGSHHHHHH

More info

UNIPROT ID

UNIPROT LINK

SYNONYMS

For research use only.