

## Human hCSE Recombinant Protein

PX-P4065-10

### OVERVIEW

---

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	Full-length
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	46.56kDa
DELIVERY CONDITION	

### PRODUCT INFORMATION

---

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	0,75
PROTEIN ACCESSION	
FORM	
BUFFER	
STABILITY & STORAGE	

### PROTEIN SEQUENCE

---

MGSHHHHHHSGLEVLFGQPMQEKDASSQGFLPHFQHFATQAIHVGQDPEQWTSRAVVPISLS  
TTFKQGAPGQHSQFEYSRSGNPTRNCLEKAVAALDGAKEYCLAFASGLAATVTITHLLKAGDQI  
ICMDDVYGGTNRVFRQVASEFGLKISFVDCSKIKLLEAAITPETKLVWIETPTNPTQKVIDIEGC  
AHIVHKHGDIIIVDNTFMSPYFQRPLALGADISMYSATKYMNGHSDVVMGLVSVNCELSHN  
RLRFLQNSLGAVPSPIDCYLCNRGLKTLHVRMEKHFKNGMAVAQFLESNPWVEKVIYPGLPS  
HPQHELVKRQCTGCTGMVTFYIKGTLQHAEIFLKNLKLFTLAESLGGFESLAELPAIMTHASVL  
KNDRDVLGISDTLIRLSVGLLEDEEDLLEDLDQALKAHPPSGSHS

### More info

UNIPROT ID

UNIPROT LINK

SYNONYMS

**For research use only.**