

Human hTEAD1 Recombinant Protein

PX-P4066-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	Partial
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	51.68kDa
DELIVERY CONDITION	

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	0,90000000000000002
PROTEIN ACCESSION	
FORM	
BUFFER	
STABILITY & STORAGE	

PROTEIN SEQUENCE

GSRSIGTTKLRLVEFSAFLEQQRDPDSYNKHLFVHIGHANHSYSDP LLESVDIRQIYDKFPEKKG
GLKELFGKGPQNAFFLVKFWADLNCNIQDDAGAFYGVTSQYESSENMTVTCSTKVCSFGKQV
VEKVETEARFENGRFVYRINRSPMCEYMINFIHKLKHLPEKYMMNSVLENFTILLVVTNRDT
QETLLCMACVFEVSNSEHGAQH HIYRLVKD

More info

UNIPROT ID

UNIPROT LINK

SYNONYMS

For research use only.