

Clostridium chauvoei CCTA-His Recombinant Protein

PX-P4067-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Clostridium chauvoei
FRAGMENT	Partial
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	33.9kDa
DELIVERY CONDITION	

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	0,900000000000000002
PROTEIN ACCESSION	
FORM	
BUFFER	
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-47% glycerol improves cryoprotection)

PROTEIN SEQUENCE

MNHKVVHHHHHHMQENTCIVETPSEGVKTFSTSDTAYADYNCFKTNLSVTFIEDQHNNQLT
ALVSTEGSFIPSGLSRVGGYYQADMYWPSKYTTLLTTYDRNNRVKITKSIPTNQIDTVSVSE
TMGYSIGGSLSIEYGKEGPKAGGGINGSYTAQRSVTYDQPDYRLLMKDSVNSASWEVAF
NATKDGYDRDSYHGIYGNQLFMRYRLYNTGINNLTDDNNLSSLIVGGFSPKVVIALTAPKG
TEESTVKVEYNRFNDQYRLRWSGTEWYGENNRNSRIDSSSESFILNWKNHHTVEHAGY

For research use only.

More info

UNIPROT ID

I6VBE8

UNIPROT LINK

<http://www.uniprot.org/uni/prot/I6VBE8>

SYNONYMS