

NME1-non-metastatic cells 1, protein (NM23A) expressed in Gene

Reference: PTXBC000293

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AWD; GAAD; NBS; NDKA; NDPK-A; NDPKA; NM23; NM23-H1; nucleoside diphosphate kinase A; NDP kinase A; granzyme A-activated DNase; metastasis inhibition factor nm23; non-metastatic cells 1, protein (NM23A) expressed in; tumor metastatic process-associated pro		

For research use only.

GENE DESCRIPTION	non-metastatic cells 1, protein (NM23A) expressed in
GENE ID	4830
ORIGIN SPECIES	Human
ACCESSIONS	BC000293
VECTOR	pDONR223
NCBI SYMBOL	NME1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaactgtgagcgtacattcattgcgatcaaacagatgggggtccagcgggggtctgtgggagagattatcaagcgtttgagcagaaaggattccgcttgggtctgaaattcatgcaagctccgaagatcttcaaggaacactacgttgacctgaaggaccgtccattctttgc
cggcctgggtgaatacatgcactcagggccggtagttgccatgggtctgggaggggctgaatgtgtggaagacgggcccagtgatgctcggggagaccaacccctgcagactccaagcctgggaccatccgtggagactctcgcatacaagttggcaggaacattatcatggcagtgatt
ctgtggagagtgacagagaagagatcggcttgggtttcacctgaggaactggtgattacacgagctgtgctcagaactggatcatgaatga