

TEX261-testis expressed 261 Gene

Reference: PTXBC020251

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein TEX261; TEG-261; testis expressed sequence 261; testis expressed 261		
For research use only.			

GENE DESCRIPTION	testis expressed 261
GENE ID	113419
ORIGIN SPECIES	Human
ACCESSIONS	BC020251
VECTOR	pDONR223
NCBI SYMBOL	TEX261
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatctggttccaccgctgtactgattggcctctacgtcttggagcgttccccaccagcatgattggagtggtgctattcaccacccctctctactttggcctctccagaccttccccctcatcatgctgacctgcctaaacttcatctgtgtggtgactagtgggtggaatcattacctag
catttcagtttttgcagaagaataattatccctctcagaggtctctggcctatttctctctgtgtgataattcgttttgcgtttttgtgtcactttcgccgggggagaacgtctgccccttaccatgcagccaggagatgatgtcttccaattattcaccaaaggcaagcggggcaaacg
cttagggatcctggtgtgtctctctcctcaaaaggccattctaccagtcgtcagaagataactga