

Human Interleukin-1 receptor type 1 recombinant protein

PX-P4006-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	64.3kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,9499999999999996
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

MKVLRLRLICFIALLISSLEADKCKEREKIKILVSSANEIDVRPCPLNPNHKGTTITWYKDDSKTPVSTEQASRIHQHKEKLVFV
PAKVEDSGHYCVVRNSSYCLRIKISAKFVENEPNLCYNAQAIFKQKLPVAGDGGGLVCPYMEFFKNENNELPKLQWYKDC
KPLLLDNIHFSGVKDRLIVMNVAEKHRGNYTCHASYTYLGKQYPITRVIEFITLEENKPTRPVIVSPANETMEVDLGSQIQLI
CNVTGQLSDIAYWKWNGSVIDEDPPVLGEDIYSVENPANKRRSTLITVLNISEIESRFYKHPFTCFKANTHGIDAAYIQLIYP
VTDDDDKDKTHTCPCPAPELLGGPSVFLFPPKPKDQLMISRTPVETCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPRE
EQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTKAKGQPREPQVYVTLPPSREEMTKNQVSLTCLVKGFYP
SDIAVEWESNGQPENNYKTTTPVLDSDSGSFFLYSKLTVDKSRWQQGNVFCSSVLHEALHNHYTQKSLSLSPGK

More info

UNIPROT ID

P14778

UNIPROT LINK

<http://www.uniprot.org/uni-prot/P14778>

SYNONYMS

IL1R1, Interleukin-1 receptor type 1, IL-1R-1, CD121a, CD121-A, IL1-R1, IL1R, IL1RA, P80, CD121 Antigen-Like Family Member A, Interleukin-1 Receptor Alpha