

Human CD278/ICOS recombinant protein

PX-P4035-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	42.34kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,6999999999999996
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	50 mM Tris-HCl pH 8, 150 mM NaCl 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MKSGLWYFFLFLCLRIKVLVTGEINGSANYEMFIFHNGGVQILCKYPDIVQQFKMQLLKGGQILC
DLTKTKGSGNTVSIKSLKFCHSQLSNNSVSFFLYNLDHSHANYFCNLSIFDPPPFPKVTLTGGY
LHIYESQLCCQLKSRDDDDKDKTHTCPPCPAPPELLGGPSVFLFPPKPKDQLMISRTPEVTCVTV
DVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVVS
NKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPE
NNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFCSSLVHEALHNHYTQKSLSLSPGK

More info

UNIPROT ID

Q9Y6W8

UNIPROT LINK

[http://www.uniprot.org/uni
prot/Q9Y6W8](http://www.uniprot.org/uni
prot/Q9Y6W8)

SYNONYMS

CD278, AILIM, Inducible
T-Cell Costimulator,
Activation-Inducible
Lymphocyte
Immunomediatory
Molecule